

## Carbohydrate Sulfotransferase 3

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

|               |   |
|---------------|---|
| Description   | Carbohydrate Sulfotransferase 3 (EC:2.8.2.17) catalyzes the transfer of sulfate groups to the 6-position of GalNAc residues in chondroitin. |
| Targets(IC50) | Endogenous Metabolite   |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481