


α -D-Glucopyranosyl β -D-arabino-2-hexulofuranosidonic acid

Chemical Properties

CAS No. : 150787-99-8
 Formula: C₁₂H₂₀O₁₂
 Molecular Weight: 356.28
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	α -D-Glucopyranosyl β -D-arabino-2-hexulofuranosidonic acid is a reagent used in biochemical reactions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8068 mL	14.0339 mL	28.0678 mL
5 mM	0.5614 mL	2.8068 mL	5.6136 mL
10 mM	0.2807 mL	1.4034 mL	2.8068 mL
50 mM	0.0561 mL	0.2807 mL	0.5614 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481