


ImmTher

Chemical Properties

CAS No. :	143005-30-5
Formula:	C ₆₅ H ₁₁₇ N ₇ O ₂₀
Molecular Weight:	1316.66
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	ImmTher is a liposome-encapsulated lipophilic disaccharide tripeptide derivative of muramyl dipeptide. It can upregulate specific monocyte cytokine genes and activate monocyte-mediated tumoricidal activity.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7595 mL	3.7975 mL	7.595 mL
5 mM	0.1519 mL	0.7595 mL	1.519 mL
10 mM	0.0759 mL	0.3797 mL	0.7595 mL
50 mM	0.0152 mL	0.0759 mL	0.1519 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481