

2,3,5-Tri-O-benzyl- β -L-arabinofuranose

Chemical Properties


CAS No. : 89615-42-9

Formula: C₂₆H₂₈O₅

Molecular Weight: 420.5

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	2,3,5-Tri-O-benzyl- β -L-arabinofuranose is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3781 mL	11.8906 mL	23.7812 mL
5 mM	0.4756 mL	2.3781 mL	4.7562 mL
10 mM	0.2378 mL	1.1891 mL	2.3781 mL
50 mM	0.0476 mL	0.2378 mL	0.4756 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481