


## trans-2-Tridecen-1-ol

## Chemical Properties

CAS No. : 68480-25-1  
 Formula: C<sub>13</sub>H<sub>26</sub>O  
 Molecular Weight: 198.35  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	trans-2-Tridecen-1-ol is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.0416 mL	25.208 mL	50.4159 mL
5 mM	1.0083 mL	5.0416 mL	10.0832 mL
10 mM	0.5042 mL	2.5208 mL	5.0416 mL
50 mM	0.1008 mL	0.5042 mL	1.0083 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481