


## Ethenyl 10-undecenoate

## Chemical Properties

CAS No. : 5299-57-0  
 Formula: C<sub>13</sub>H<sub>22</sub>O<sub>2</sub>  
 Molecular Weight: 210.31  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Ethenyl 10-undecenoate is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.7549 mL	23.7744 mL	47.5489 mL
5 mM	0.951 mL	4.7549 mL	9.5098 mL
10 mM	0.4755 mL	2.3774 mL	4.7549 mL
50 mM	0.0951 mL	0.4755 mL	0.951 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481