

4-Decenoic Acid

Chemical Properties


CAS No. : 26303-90-2

Formula: C₁₀H₁₈O₂

Molecular Weight: 170.25

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	4-Decenoic Acid is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.8737 mL	29.3686 mL	58.7372 mL
5 mM	1.1747 mL	5.8737 mL	11.7474 mL
10 mM	0.5874 mL	2.9369 mL	5.8737 mL
50 mM	0.1175 mL	0.5874 mL	1.1747 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481