

cis-3-Octen-1-ol

## Chemical Properties

CAS No. : 20125-84-2  
 Formula: C<sub>8</sub>H<sub>16</sub>O  
 Molecular Weight: 128.21  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	cis-3-Octen-1-ol is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.7997 mL	38.9985 mL	77.997 mL
5 mM	1.5599 mL	7.7997 mL	15.5994 mL
10 mM	0.780 mL	3.8999 mL	7.7997 mL
50 mM	0.156 mL	0.780 mL	1.5599 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481