


Methyl lesquerolate

Chemical Properties

CAS No. :	4102-96-9
Formula:	C ₂₁ H ₄₀ O ₃
Molecular Weight:	340.54
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Methyl lesquerolate is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9365 mL	14.6826 mL	29.3651 mL
5 mM	0.5873 mL	2.9365 mL	5.873 mL
10 mM	0.2937 mL	1.4683 mL	2.9365 mL
50 mM	0.0587 mL	0.2937 mL	0.5873 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481