

2,3,4,6-Tetra-O-benzyl-D-mannopyranose


Chemical Properties

CAS No. : 61330-61-8

Formula: C₃₄H₃₆O₆

Molecular Weight: 540.65

Storage: Store under nitrogen
 Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

| | |
|---------------|---|
| Description | 2,3,4,6-Tetra-O-benzyl-D-mannopyranose has potential applications in glycobiology research. |
| Targets(IC50) | Others |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|------------|
| 1 mM | 1.8496 mL | 9.2481 mL | 18.4963 mL |
| 5 mM | 0.3699 mL | 1.8496 mL | 3.6993 mL |
| 10 mM | 0.185 mL | 0.9248 mL | 1.8496 mL |
| 50 mM | 0.037 mL | 0.185 mL | 0.3699 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481