


p-Nitrophenyl 2-acetamido-2-deoxy-3,6-di-O-acetyl-β-D-glucopyranoside

Chemical Properties

CAS No. : 61891-87-0
 Formula: C₁₈H₂₂N₂O₁₀
 Molecular Weight: 426.38
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	p-Nitrophenyl 2-acetamido-2-deoxy-3,6-di-O-acetyl-β-D-glucopyranoside is a reagent used in biochemical reactions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3453 mL	11.7266 mL	23.4533 mL
5 mM	0.4691 mL	2.3453 mL	4.6907 mL
10 mM	0.2345 mL	1.1727 mL	2.3453 mL
50 mM	0.0469 mL	0.2345 mL	0.4691 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481