


Sulfo-SANPAH sodium

Chemical Properties

CAS No. :	144672-60-6
Formula:	C ₁₆ H ₁₇ N ₆ NaO ₉ S
Molecular Weight:	492.4
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Sulfo-SANPAH sodium is a reagent used in biochemical reactions.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0309 mL	10.1543 mL	20.3087 mL
5 mM	0.4062 mL	2.0309 mL	4.0617 mL
10 mM	0.2031 mL	1.0154 mL	2.0309 mL
50 mM	0.0406 mL	0.2031 mL	0.4062 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481