


(R)-3-Hydroxyisobutyric acid sodium

Chemical Properties

CAS No. :	1228078-57-6
Formula:	C ₄ H ₇ NaO ₃
Molecular Weight:	126.09
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(R)-3-Hydroxyisobutyric acid sodium is a biochemical marker of the metabolic disorder, 3-hydroxyisobutyric aciduria, characterized by increased concentrations in urine.
Targets(IC50)	Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.9308 mL	39.6542 mL	79.3084 mL
5 mM	1.5862 mL	7.9308 mL	15.8617 mL
10 mM	0.7931 mL	3.9654 mL	7.9308 mL
50 mM	0.1586 mL	0.7931 mL	1.5862 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481