


(Rac)-1,2-Dihexadecylglycerol

Chemical Properties

CAS No. : 13071-60-8
 Formula: C₃₅H₇₂O₃
 Molecular Weight: 540.94
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	(Rac)-1,2-Dihexadecylglycerol is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8486 mL	9.2432 mL	18.4863 mL
5 mM	0.3697 mL	1.8486 mL	3.6973 mL
10 mM	0.1849 mL	0.9243 mL	1.8486 mL
50 mM	0.037 mL	0.1849 mL	0.3697 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481