


1,4-Cyclohexanedimethanol Dilaurate

Chemical Properties

CAS No. : 135025-34-2
 Formula: C₃₂H₆₀O₄
 Molecular Weight: 508.82
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	1,4-Cyclohexanedimethanol Dilaurate is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9653 mL	9.8267 mL	19.6533 mL
5 mM	0.3931 mL	1.9653 mL	3.9307 mL
10 mM	0.1965 mL	0.9827 mL	1.9653 mL
50 mM	0.0393 mL	0.1965 mL	0.3931 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481