


Benzyl 2-deoxy-2-phthalimido-3,6-di-O-benzyl-β-D-glucopyranoside

Chemical Properties

CAS No. : 80035-36-5
 Formula: C₃₅H₃₃N₁O₇
 Molecular Weight: 579.64
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Benzyl 2-deoxy-2-phthalimido-3,6-di-O-benzyl-β-D-glucopyranoside is a reagent used in biochemical reactions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7252 mL	8.626 mL	17.2521 mL
5 mM	0.345 mL	1.7252 mL	3.4504 mL
10 mM	0.1725 mL	0.8626 mL	1.7252 mL
50 mM	0.0345 mL	0.1725 mL	0.345 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481