


## Poly(ethylene glycol) Monooleate (MW 860)

### Chemical Properties

CAS No. : 9004-96-0  
 Formula: C<sub>20</sub>H<sub>38</sub>O<sub>3</sub>  
 Molecular Weight: 326.51  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



### Biological Description

Description	Poly(ethylene glycol) Monooleate (MW=860) functions as a plasticizer and surfactant.
Targets(IC50)	Others

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0627 mL	15.3135 mL	30.6269 mL
5 mM	0.6125 mL	3.0627 mL	6.1254 mL
10 mM	0.3063 mL	1.5313 mL	3.0627 mL
50 mM	0.0613 mL	0.3063 mL	0.6125 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481