


5-Bromo-4-chloro-3-indolyl- α -L-fucopyranoside

Chemical Properties

CAS No. :	171869-92-4
Formula:	C ₁₄ H ₁₅ BrClNO ₅
Molecular Weight:	392.63
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-Bromo-4-chloro-3-indolyl- α -L-fucopyranoside can serve as a chromogenic substrate for α -L-Fucosidase.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5469 mL	12.7346 mL	25.4693 mL
5 mM	0.5094 mL	2.5469 mL	5.0939 mL
10 mM	0.2547 mL	1.2735 mL	2.5469 mL
50 mM	0.0509 mL	0.2547 mL	0.5094 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481