

Oleyl oleate

Chemical Properties


CAS No. : 3687-45-4

Formula: C₃₆H₆₈O₂

Molecular Weight: 532.92

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Oleyl oleate is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8765 mL	9.3823 mL	18.7645 mL
5 mM	0.3753 mL	1.8765 mL	3.7529 mL
10 mM	0.1876 mL	0.9382 mL	1.8765 mL
50 mM	0.0375 mL	0.1876 mL	0.3753 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481