


## S-(3)-Hydroxymyristic acid, methyl ester

## Chemical Properties

CAS No. :	76835-67-1
Formula:	C15H30O3
Molecular Weight:	258.4
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	S-(3)-Hydroxymyristic acid, methyl ester is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.870 mL	19.3498 mL	38.6997 mL
5 mM	0.774 mL	3.870 mL	7.7399 mL
10 mM	0.387 mL	1.935 mL	3.870 mL
50 mM	0.0774 mL	0.387 mL	0.774 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481