


18:1 Hemi BMP (S,R)

Chemical Properties

CAS No. :	474943-36-7
Formula:	C60H114NO11P
Molecular Weight:	1056.52
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	18:1 Hemi BMP (S,R) is a reagent used in biochemical reactions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9465 mL	4.7325 mL	9.465 mL
5 mM	0.1893 mL	0.9465 mL	1.893 mL
10 mM	0.0947 mL	0.4733 mL	0.9465 mL
50 mM	0.0189 mL	0.0947 mL	0.1893 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481