

DSPE-PEG-FITC (MW 2000)

Chemical Properties

CAS No. : 1000788-53-3

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG-FITC, MW 2000 is a fluorescently labeled PEG lipid designed for liposome preparation as a drug delivery vehicle for targeted therapy. FITC, a green dye, has a peak absorption at 494 nm and maximum emission at 520 nm, making it suitable for staining biological samples or nanoparticles.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481