

Ammonium polyphosphate

Chemical Properties

CAS No. : 68333-79-9
Formula: $((\text{NH}_4)_3\text{PO}_4)_n$
Molecular Weight:
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Ammonium polyphosphate, n < 20, water solubility > 90g/100ml is a reagent used in biochemical reactions.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481