


Phytic acid calcium

Chemical Properties

CAS No. :	23183-60-0
Formula:	C ₆ H ₁₈ CaO ₂₄ P ₆
Molecular Weight:	700.11
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Phytic acid calcium is a reagent used in biochemical reactions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4283 mL	7.1417 mL	14.2835 mL
5 mM	0.2857 mL	1.4283 mL	2.8567 mL
10 mM	0.1428 mL	0.7142 mL	1.4283 mL
50 mM	0.0286 mL	0.1428 mL	0.2857 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481