

TMA chloride, for molecular biology

Chemical Properties


CAS No. : 75-57-0

Formula: C₄H₁₂ClN

Molecular Weight: 109.6

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	TMA chloride, for molecular biology, 99%(T) (Tetramethylammonium chloride, for molecular biology, 99%(T)), enhances hybridization specificity and increases the DNA melting temperature in PCR.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	9.1241 mL	45.6204 mL	91.2409 mL
5 mM	1.8248 mL	9.1241 mL	18.2482 mL
10 mM	0.9124 mL	4.562 mL	9.1241 mL
50 mM	0.1825 mL	0.9124 mL	1.8248 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481