


Sodium 1-octanesulfonate monohydrate

Chemical Properties

CAS No. : 207596-29-0
 Formula: C₈H₁₉NaO₄S
 Molecular Weight: 234.29
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Sodium 1-octanesulfonate monohydrate, 97%(T), is an ion-pairing agent characterized by anionic surfactant properties. It interacts with protonated amino groups in Tobramycin molecules.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.2682 mL	21.3411 mL	42.6821 mL
5 mM	0.8536 mL	4.2682 mL	8.5364 mL
10 mM	0.4268 mL	2.1341 mL	4.2682 mL
50 mM	0.0854 mL	0.4268 mL	0.8536 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481