

Sodium polyphosphate

Chemical Properties

CAS No. :	68915-31-1
Formula:	(NaPO ₃) _n
Molecular Weight:	
Storage:	Store at RT <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Sodium polyphosphate, crystal, +200 mesh, 96%, is a water-soluble inorganic stabilizer that functions as both a surfactant and a metal ion chelating agent. It is widely used in food, cosmetics, and industrial products to enhance texture and stability. Additionally, it helps prevent degradation caused by metal ions, thus extending product shelf life.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481