

EDTA diammonium salt hydrate

Chemical Properties

CAS No. : 304675-80-7


Formula: C₁₀H₂₄N₄O₉

Molecular Weight: 344.32

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	EDTA diammonium salt hydrate, 97%, is a derivative of ethylenediaminetetraacetic acid (EDTA) and is widely used as a chelating agent because of its ability to bind with various metal ions to form stable complexes.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9043 mL	14.5214 mL	29.0428 mL
5 mM	0.5809 mL	2.9043 mL	5.8086 mL
10 mM	0.2904 mL	1.4521 mL	2.9043 mL
50 mM	0.0581 mL	0.2904 mL	0.5809 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481