


## Sucrose monodecanoate

## Chemical Properties

CAS No. :	31835-06-0
Formula:	C <sub>22</sub> H <sub>42</sub> O <sub>13</sub>
Molecular Weight:	514.56
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Sucrose monodecanoate is a biochemical reagent that can be used as a biomaterial for life science related research and as a sulfonylation reagent for organic synthesis and drug discovery.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9434 mL	9.717 mL	19.4341 mL
5 mM	0.3887 mL	1.9434 mL	3.8868 mL
10 mM	0.1943 mL	0.9717 mL	1.9434 mL
50 mM	0.0389 mL	0.1943 mL	0.3887 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481