

## Luviquat FC 550

## Chemical Properties

CAS No. : 95144-24-4  
Formula:  $(C_6H_9N_2.C_6H_9NO.Cl)_x$   
Molecular Weight:  
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Luviquat FC 550 is a cationic polymer frequently utilized in personal care and cosmetic products. It is a viscous liquid containing quaternary ammonium groups, giving it a positive charge. Luviquat FC 550 exhibits properties suitable for these applications, such as conditioning, anti-static, and film-forming capabilities for hair and skin. Additionally, it serves as a binder in oral care products and a flocculant in industrial water treatment.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481