


Sorbitan Sesquioleate

Chemical Properties

CAS No. : 8007-43-0
Formula: C₆₆H₁₂₆O₁₆
Molecular Weight: 1175.7
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Sorbitan Sesquioleate is a synthetic surfactant and emulsifier derived from sorbitol and modified oleic acid.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8506 mL	4.2528 mL	8.5056 mL
5 mM	0.1701 mL	0.8506 mL	1.7011 mL
10 mM	0.0851 mL	0.4253 mL	0.8506 mL
50 mM	0.017 mL	0.0851 mL	0.1701 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481