

ι-Carrageenan

Chemical Properties

CAS No. : 9062-07-1

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	ι-Carrageenan (Viscarin SD 309) is a biochemical reagent. It can be isolated from certain red algae species, specifically from the genera Eucheuma and Chondrus. It has potential applications in the flocculation and stabilization of protein emulsions.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481