


Fast Red RC salt

Chemical Properties

CAS No. : 85252-22-8
 Formula: C₁₄H₁₂Cl₆N₄O₂Zn
 Molecular Weight: 546.4
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Fast Red RC salt, 97%, is utilized in plant science to detect aminotransferase activity and serves as a marker for phosphatase activity in diagnostic procedures.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8302 mL	9.1508 mL	18.3016 mL
5 mM	0.366 mL	1.8302 mL	3.6603 mL
10 mM	0.183 mL	0.9151 mL	1.8302 mL
50 mM	0.0366 mL	0.183 mL	0.366 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481