


D,L-Mevalonic acid dicyclohexylammonium

Chemical Properties

CAS No. : 1215802-31-5
 Formula: C₁₈H₃₅NO₄
 Molecular Weight: 329.48
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	D,L-Mevalonic acid dicyclohexylammonium is a reagent used in biochemical reactions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0351 mL	15.1754 mL	30.3509 mL
5 mM	0.607 mL	3.0351 mL	6.0702 mL
10 mM	0.3035 mL	1.5175 mL	3.0351 mL
50 mM	0.0607 mL	0.3035 mL	0.607 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481