


## 7-Epi lincomycin hydrochloride

## Chemical Properties

CAS No. :	26389-84-4
Formula:	C <sub>18</sub> H <sub>35</sub> ClN <sub>2</sub> O <sub>6</sub> S
Molecular Weight:	443
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	7-Epi lincomycin hydrochloride is a reagent used in biochemical reactions.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2573 mL	11.2867 mL	22.5734 mL
5 mM	0.4515 mL	2.2573 mL	4.5147 mL
10 mM	0.2257 mL	1.1287 mL	2.2573 mL
50 mM	0.0451 mL	0.2257 mL	0.4515 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481