

m-PEG-Vinylsulfone (MW 20000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	m-PEG-Vinylsulfone (MW 20000) (m-PEG-VS (MW 20000)) is a PEG derivative that modulates drug release in parenteral delivery systems. It is applicable in drug delivery research.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481