

PLLA4000-PEG6000-PLLA4000

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	PLLA4000-PEG6000-PLLA4000 is an amphiphilic triblock polymer derived from polylactic acid (PLA), designed to enhance the specificity and cell affinity of PLA-based biomaterials. This compound can be utilized in studies related to drug delivery.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481