

PLLA10000-PEG2000-RhB

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PLLA10000-PEG2000-RhB is a polylactic acid derivative with a terminal rhodamine label (PLLA10000-PEG2000-RhB). It is used for labeling nanomicelles to enable real-time monitoring of drug release within the body.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481