

## Poly(ethylene glycol) methyl ether methacrylate (MW 950)

### Chemical Properties

CAS No. : 26915-72-0

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Poly(ethylene glycol) methyl ether methacrylate (MW 950) is a non-linear analog of polyethylene glycol (PEG) and can be utilized for PEG surface modification.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481