

DOPE-PEG-COOH (MW 1000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DOPE-PEG-COOH (MW 1000) is a phospholipid characterized by excellent biocompatibility and notable amphiphilic properties, making it an ideal and pivotal excipient or dosage form in pharmaceutical formulations for enhanced therapeutic outcomes. It is applicable in studies related to drug delivery.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481