

DOPE-PEG-COOH (MW 3400)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DOPE-PEG-COOH (MW 3400) is a phospholipid known for excellent biocompatibility and significant amphiphilic properties, making it an ideal and essential excipient or form in drug formulations for enhanced therapeutic efficacy. It is applicable in the research of drug delivery.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481