

DSPE-PEG-OH (MW 3400)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DSPE-PEG-OH, with a molecular weight of 3400, is an amphiphilic and polydisperse PEG that spontaneously self-assembles in aqueous solutions. This polymer is utilized in the formulation of liposomes for targeted drug delivery. Reagent grade, for research use only.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481