

DSPE-PEG-SH (MW 3400)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG-SH, MW 3400, is a maleimide-reactive PEG lipid. This amphiphilic polymer can form lipid bilayers in water and is utilized in the creation of liposomes. It serves as a drug nanocarrier for delivering mRNA, DNA, or siRNA. Reagent grade, intended for research purposes only.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481