

DSPE-PEG-IA (MW 5000)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DSPE-PEG-DBCO, MW 2000, is a phospholipid PEG polymer that incorporates cyclooctyne. This polymer spontaneously self-assembles into micelles or lipid bilayers in aqueous solutions, making it suitable for nanoparticle or liposome preparation as drug carriers in targeted drug delivery systems. DBCO can react with azide molecules via copper-free click chemistry to form stable triazole bonds. Reagent-grade, for research use only.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481