

DSPE-PEG-Cy3, MW 2000

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG-Cy3, MW 2000 is a PEG lipid conjugated with a fluorescent group. The Cy3 fluorophore is commonly utilized in applications such as immunolabeling, nucleic acid labeling, fluorescence microscopy, and flow cytometry. This dye exhibits peak absorption at wavelengths between 548-552 nm and maximum emission between 562-570 nm.
Targets(IC50)	Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481