

DSPE-PEG-Cy5, MW 2000

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Keep away from moisture,Keep away from direct sunlight

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG-Cy5, MW 2000, is a PEG lipid conjugated with a fluorescent group. The Cy5 fluorophore is commonly applied in the labeling of proteins and nucleic acids, suitable for imaging, flow cytometry, and genomic applications. The dye features an excitation peak at 651 nm and an emission peak at 670 nm.
Targets(IC50)	Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481