

DSPE-PEG-Cy5, MW 3400

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG-Cy5, MW3400, is a PEG lipid conjugated with a fluorescent group. The Cy5 fluorophore is commonly utilized for labeling proteins and nucleic acids, facilitating imaging, flow cytometry, and genomic applications. The dye exhibits an excitation peak at 651 nm and an emission peak at 670 nm.
Targets(IC50)	Liposome

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481