

DSPE-PEG-Cy5, MW 5000

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DSPE-PEG-Cy5, MW5000, is a PEG lipid conjugated with a fluorescent group. The Cy5 fluorophore is commonly utilized for protein and nucleic acid labeling for imaging, flow cytometry, and genomic applications. The dye exhibits an excitation peak at 651 nm and an emission peak at 670 nm.
Targets(IC50)	Liposome

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481