

## Stearic acid-PEG-Rhodamine (MW 2000)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Stearic acid-PEG-Rhodamine, MW 2000, is an amphiphilic PEG lipid capable of forming micelles in water. Rhodamine provides staining capabilities for samples and facilitates tracking with fluorescence microscopy, with a maximum absorption wavelength of 570 nm and an emission peak at 595 nm.
Targets(IC50)	Liposome

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481