

Stearic acid-PEG-Mal (MW 1000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Stearic acid-PEG-Mal (MW 1000) is a click reagent that can form triazole bonds by reacting with terminal alkynes through CuAAC or with cyclooctyne (DBCO/BCN) via SPAAC. Silane.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481